

<b>Description</b>	<b>Precipitating polyclonal swine antiserum to rabbit light chains, free and bound</b>	
<b>Product code</b>	SwAR/L(SD+HD)	
<b>Biological origin</b>	Swine	
<b>Physical form</b>	Delipidated, heat inactivated lyophilized stable whole serum	
<b>Preservative</b>	No preservative added, as it may interfere with the antibody activity. No foreign protein added.	
<b>Immunogen</b>	Purified light chains isolated from polyclonal rabbit IgG. Freund's complete adjuvant is used in the first step of the immunization procedure.	
<b>Adsorption</b>	Immunoaffinity adsorbed using insolubilized antigens as required to eliminate antibodies reacting with the Fc fragments of immunoglobulins or with other plasma proteins.	
<b>Identity &amp; Specificity</b>	As tested in immunoelectrophoresis and double radial immunodiffusion against pooled normal rabbit serum and purified IgG fractions and light chain preparations, the antiserum reacts with determinants exposed on the surface (SD) and the Fab portion of all intact immunoglobulins, with Fab or F(ab') <sub>2</sub> subunits, and with free light chains of both types; also with determinants exposed only on free light chains but hidden in the intact immunoglobulin molecule (HD or hidden determinants) or Bence Jones proteins. This antiserum is not kappa or lambda chain specific. It does not react with any non-Ig protein in rabbit serum, as tested by immunoelectrophoresis and double radial immunodiffusion.	
<b>Cross-reactivity</b>	Inter-species cross-reactivity is a normal feature of antibodies to immunoglobulins, since Ig of different species frequently share antigenic determinants. Cross-reactivity of this antiserum has not been tested in detail, but a wide scale of reactivity with IgG of other mammalian species must be assumed. The degree of cross-reactivity depends in addition to the phylogenetic relationship also on the concentrations of the reactants and the sensitivity of the assay format.	
<b>Protein concentration</b>	Total protein and IgG concentration in the antiserum are comparable to those of pooled swine serum. No foreign proteins added.	
<b>Antibody titre</b>	Precipitin titre not less than 1:64 when tested against normal rabbit serum in agar block titration.	
<b>Intended use</b>	As precipitating antiserum in a variety of immunodiffusion techniques including immunoelectrophoresis and double and single radial immunodiffusion (Mancini, Ouchterlony). It has not been tested for use in nephelometry, ELISA or immunochemistry, but this does not exclude such use if proper controls are included.	
<b>Directions for use</b>	In immunoelectrophoresis use 2 µl or equivalent against 120 µl antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 µl antiserum in a 3 mm diameter centre well and 2 µl serum samples (neat and serially diluted) in 2 mm diameter peripheral wells. In single immunodiffusion, use 1% antiserum in agar gel to establish a standard curve.	
<b>Handling</b>	The lyophilized antiserum is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. It is reconstituted by adding 1 ml sterile distilled water. Working dilutions are prepared by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product.	
<b>Packing</b>	Vial with 1 ml lyophilized antiserum.	
<b>Storage / shelf life</b>	Lyophilized at +4°C	at least 10 years
	reconstituted at or below -20°C	3-5 years
	reconstituted at +4°C	7 days
<b>Caution</b>	This antiserum should be handled by qualified persons only and appropriate precautions should be taken in its handling and disposal, and of all associated materials. For <i>in vitro</i> laboratory research purposes only.	

**NORDIC IMMUNOLOGICAL LABORATORIES**  
Langendijk 25, 5652 AX Eindhoven, The Netherlands  
Tel. +31 630 070 625, Fax: +31 402 920 069  
E-mail: [info@nordiclabs.nl](mailto:info@nordiclabs.nl)  
[www.nordiclabs.nl](http://www.nordiclabs.nl)